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AUTOMATIC PARTICLE DETECTION AND LEAK TEST EQUIPMENT

FACTORY ACCEPTANCE TEST

Machine **Automatic Inspection and Leak Test Machine**

Model **RS1**

Serial No. **618**

Customer **TFDA**

FAT RS1	TFDA	Serial No: 618	Pag. 1 - 17
---------	------	----------------	-------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 -
C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133
website : www.cmp-pharma.it email : info@cmp-pharma.it

1. EQUIPMENT DETAILS

Equipment Type: Automatic Inspection and Leak Test Machine
Model: RS1
Manufacturer: CMP PHAR.MA S.R.L.
Manufacturer Details: Via G. Galilei, 26
36030 Costabissara (Vicenza) - Italy
TEL: +39 (0444) 971716
FAX: +39 (0444) 553262

2. FAT DETAILS

Participants of customer:

張 翕 偉
William Cheng
Jason Chen

Participants of CMP PHAR.MA S.R.L.:

Daniela Tassou

Date of F.A.T.:

2017.08.28

FAT RS1	TFDA	Serial No: 618	Pag. 2 - 17
---------	------	----------------	-------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 -
C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133
website : www.cmp-pharma.it email : info@cmp-pharma.it

GENERAL:

3.1 Purpose

The purpose of this document is to describe the FAT (Factory Acceptance Test) that will be followed by the Customer on the fully automatic inspection and leak test machine mod. RS1 manufactured by CMP PHAR.MA S.R.L.

3.2 References

1. CMP PHAR.MA S.R.L. Instruction User's Manual of the machine

3.3 Responsibility

CMP PHAR.MA S.R.L. representatives are responsible for the review of this document, they will verify the records and write the report.
The Customer personnel assigned to conduct the FAT will oversee the process, perform the qualification and record the information.
The Customer Quality Assurance is responsible for the review and approval of this protocol prior to FAT and the final FAT report.

3.4 The factory acceptance test comprises of the following:

1. Checklist of the scope of supply according to contract
2. Checklist of layout drawings / documents
3. Machine built, electrical supply, safety features, spare parts etc
4. Test runs
5. Resume of the possible modifications required
6. Dates and delivery schedule

3.5 Procedure

The tasks identified in this protocol shall be carried out as follows:

- Each task identified shall be performed.
- Each task shall be carried out as described.
- The FAT assigned personnel will decide if the acceptance criteria have been met.
- All tasks must declare either a pass or fail verdict.
- Each result/outcome must be signed and dated by the individual carrying out the test as well as a witness, if possible.

FAT RS1	TFDA	Serial No: 618	Pag. 3 - 17
---------	------	----------------	-------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 -
C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133
website : www.cmp-pharma.it email : info@cmp-pharma.it

The procedure for a failed task:

- The task procedure should be repeated.
- The task failure must be clearly noted in the comments section of the result sheet.
- Improvements or required changes must be recorded, agreed upon .

3.6 Acceptance criteria

The factory acceptance testing is complete when the test protocol has been successfully completed and the report is approved.

3.7 Completing and correcting documentation

Documents should be completed in blue ink.

All corrections must be initialed and dated with the error neatly crossed out with a single line.

FAT RS1	TFDA	Serial No: 618	Pag. 4 - 17
---------	------	----------------	-------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 –
C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA – SWIFT BPVI IT 21040 – IBAN IT 91Y0572860210040570004133
website : www.cmp-pharma.it email : info@cmp-pharma.it

4.0 CHECKLIST

4.1 Check list of the scope of supply according to contract

4.1.1 Inspection Machine model RS1

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
No. 1 camera for tip/cap control	Check station with individual camera present	Camera fitted	YES	偉	2017-08-28
No. 2 cameras for particles detection	Check stations with individual cameras present	Cameras fitted	YES	偉	2017-08-28
No. 1 camera for particle/level detection	Check stations with individual cameras present	Camera fitted	YES	偉	2017-08-28
No. 4 electrodes for leak test	Check that all parts are present	All parts are present	YES	偉	2017-08-28
Change parts suitable for ampoules: 1 ml and 2 ml	Check that all parts are present	All parts are present	YES	偉	2017-08-28

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 5 - 17
---------	------	----------------	-------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 - C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133 website : www.cmp-pharma.it email : info@cmp-pharma.it

4.1.2 Documents Assessment

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
User Manual	User Manual in English and correct version	Ensure that documents are in English and complete	YES	偉	2017.8.28
Electric and pneumatic diagram	Available in English and correct version.	Ensure that documents are in English and complete	YES	偉	2017.8.28
CE-Compliance Declaration	Available in English and correct version	Ensure that CE Declaration is in English and complete	YES	偉	2017.8.28
Spare parts list	Available in English and correct version Indicates part number and supplier	Ensure that document is in English and complete	YES	偉	2017.8.28
Mechanical and Electrical components list	Available in English and correct version Indicates part number and supplier	Ensure that document is in English and complete	YES	偉	2017.8.28
Alarms list	Available in English and correct version	Ensure that document is in English and complete	YES	偉	2017.8.28
BASLER Cameras User's Manual	Available in English and correct version	Ensure documents are in English and complete	YES	偉	2017.8.28

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 6 - 17
---------	------	----------------	-------------

4.1.3 Machine Built Assessment

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Layout plan available	Check dimensions of machine	Dimensions of the machine comply with layout plan	YES	偉	2017-08-28
Check 4 main units of equipment available	Loading unit; Inspection unit; Unloading unit and control panel unit are present	4 main units of the machine are present	YES	偉	2017-08-28
Machine is monoblock - 1 part is the machine and the other is the PC control panel	Check that the machine is monoblock	Machine is monoblock	YES	偉	2017-08-28
Starwheels for ampoules/vials feeding in POM-C plastic material	Starwheels for ampoules/vials feeding in POM-C plastic material	The ampoules/vials starwheels are in POM-C plastic material	YES	偉	2017-08-28
ORIENTAL MOTOR rotation motors	Check that No.6 installed rotation motors are ORIENTAL MOTOR (No.1 motor for slow rotation on TV1 station, No.3 motors on TV2-TV3 - TV4 stations, No.1 motor for vibrating unit, No.1 motor for leak test rotation)	No.6 installed rotation motors are ORIENTAL MOTORS	YES	偉	2017-08-28
SIEMENS motors	Check that the installed motor is SIEMENS (No.1 motor for loading conveyor belt)	The motor is SIEMENS	YES	偉	2017-08-28
SEW EURODRIVE motor	Check that main motor is SEW EURIDRIVE	The installed main motor is SEW EURODRIVE	YES	偉	2017-08-28
BASLER Cameras	Check that the N. 4 installed cameras are from BASLER	The N. 4 installed cameras are from BASLER	YES	偉	2017-08-28
Panel PC 17" touch screen	Check that the Panel PC is 17" touch screen	The panel PC is 17" touch screen	YES	偉	2017-08-28
Outfeed trays with clear segregation of good and bad products	Check that the outfeed trays have a clear segregation of good and bad products	The outfeed trays have clear segregation of good and bad products	YES	偉	2017-08-28

Comments

TFDA Ask for different colors of the labels of outfeed trays
(GREEN AND RED)

FAT RS1	TFDA	Serial No: 618	Pag. 7 - 17
---------	------	----------------	-------------

4.1.5 PLC Assessment

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
PLC (Programmable Logic Controller)	Check that PLC is working	Check run-stop led light	YES	偉	2017-03-28
	Check that some inputs are read correctly	Led lights correspond to the input state	YES	偉	2017-03-28

Comments

4.1.6 Safety Devices Assessment

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Fixed guards to stop limbs from coming into contact with moving parts	Check fixed guards present	Fixed guards present	YES	偉	2017-03-28
Protected electrical panel	Check that the electrical panel is protected	Electrical panel protected	YES	偉	2017-03-28
Hazardous sign warnings and warning signs	Check that hazardous sign warnings and warning signs are present and clearly visible	Hazardous sign warnings and warning signs are present and clearly visible	YES	偉	2017-03-28
	* Symbol warning not to remove the safety devices	Present and clearly visible			
	* Symbol warning not to work on moving parts	Present and clearly visible			
	Plate on electric panel giving voltage, phases, frequency	Present and clearly visible	YES	偉	2017-03-28

FAT RS1	TFDA	Serial No: 618	Pag. 9 - 17
---------	------	----------------	-------------

Re-settable emergency push buttons	Check that emergency push buttons are present and in accessible location	Emergency push buttons available and accessible	YES	偉	2017-01-28
Electricity isolating switch	Check that the electricity isolating switch is present and accessible	Electricity isolating switch is present and accessible	YES	偉	2017-03-28
Safety clutches	Check that safety clutches are present	Safety clutches are present	YES	偉	2017-03-28
Alarms: See the list in the technical documentation	Check that the alarms are present and operating	Alarms are present and operating	YES	偉	2017-03-28

Comments

*Labels will be placed when doors will be available

4.1.7 Emergency Stops

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Emergency stops	Press the emergency pushbuttons while the machine is running at the max. inspection speed (the test should be performed with at least 1 change part)	Machine should immediately stop	YES	偉	2017-03-28
		Possible machine part movement must not cause danger to equipment or user	YES	偉	2017-03-28
		The emergency alarm is displayed in the panel pc	YES	偉	2017-03-28
	Reset and restore start situation	The alarm is no more displayed and machine works	YES	偉	2017-03-28

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 10 - 17
---------	------	----------------	--------------

4.1.8 Software Verification

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Scada Software Version 3.0 (Supervisory control and data acquisition)	Check that the versions of the installed software is 3.0	The versions of the installed software is 3.0	YES	偉	2017.08.28
	The SCADA version should be displayed in every page	Print one screen shot and attach it to the protocol	YES	偉	2017.08.28
	The machine Serial Number should be displayed in every printed document	Print one report and attach it to the protocol	YES	偉	2017.08.28

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 11 - 17
---------	------	----------------	--------------

CMP PHARMA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 - C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133 website : www.cmp-pharma.it email : info@cmp-pharma.it

4.2 Test runs

4.2.1 Ampoules 1 ml

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Mechanical Tests are performed with the available samples	Switch on the machine, and start the test	Ampoules are conveyed smoothly through the machine	YES	偉	2017.08.29
For particle inspection, tests are completed with customer's samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for particles inspection	YES	偉	2017.08.29
Confirm maximum output of 12.000 pcs/h. <i>Note:</i> The most suitable test result to detect particles in Customer product might differ from maximum test speed.	Test machine performance at maximum speed according to the container diameter	Machine must operate efficiently during duration of testing at maximum speed	YES	偉	2017.08.29
Ampoule speed rotation can be adjusted to at least 3000 rpm	Check that ampoule speed rotation can be adjusted to at least 3000 rpm	The ampoule speed rotation can be adjusted	YES	偉	2017.08.29
For tip/cap inspection, tests are completed with customer's samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for cosmetic inspection	YES	偉	2017.08.29
For fill level inspection, tests are completed with customer's samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for fill level inspection	YES	偉	2017.08.29

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 12 - 17
---------	------	----------------	--------------

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
For leak test, certified samples provided by CMP are used	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for leak test and ampoules with certified pinholes are rejected	YES	偉	2017.08.29
Rejected ampoules are segregated into separate trays	Check that rejected ampoules are segregated into separate trays	Rejected ampoules are segregated into separate trays	YES	偉	2017.08.29
Processed ampoules are correctly counted by the machine software according to the inspection station	Check the results on machine display	Pass a known quantity of ampoules (50 or 100, for example) and check that they are counted correctly	YES	偉	2017.08.29

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 13 - 17
---------	------	----------------	--------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 - C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133 website : www.cmp-pharma.it email : info@cmp-pharma.it

4.2.2 Ampoules 2 ml

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
Mechanical Tests are performed with the available samples	Switch on the machine, and start the test	Ampoules are conveyed smoothly through the machine	YES	偉	2017.08.29
For particle inspection, tests are completed with customer samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for particles inspection	YES	偉	2017.08.29
Confirm maximum output of 12.000 pcs/h. <i>Note:</i> The most suitable test result to detect particles in Customer product might differ from maximum test speed.	Test machine performance at maximum speed according to the container diameter	Machine must operate efficiently during duration of testing at maximum speed	YES	偉	2017.08.29
For tip/cap inspection, tests are completed with customer's samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for cosmetic inspection	YES	偉	2017.08.29
For fill level inspection, tests are completed with customer's samples	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for fill level inspection	YES	偉	2017.08.29

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 14 - 17
---------	------	----------------	--------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 - C.F.P.IVA 03872190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133 website : www.cmp-pharma.it email : info@cmp-pharma.it

Prerequisites	Test Procedures	Acceptance criteria	Criteria met YES/NO	Sign	Date
For leak test, certified samples provided by CMP are used	Switch on the machine, and set the suitable inspection parameters	Set parameters are suitable for leak test	YES	偉	2017.08.29
Rejected ampoules are segregated into separate trays	Check that rejected ampoules are segregated into separate trays	Rejected ampoules are segregated into separate trays	YES	偉	2017.08.29
Processed ampoules are correctly counted by the machine software according to the inspection station	Check the results on machine display	Processed ampoules are correctly counted by the machine software according to the inspection station	YES	偉	2017.08.29

Comments

Some Certified ampoules with pinhole for leak test has been broken.
 Cmp will supply other samples.

FAT RS1	TFDA	Serial No: 618	Pag. 15 - 17
---------	------	----------------	--------------

CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 - C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133 website : www.cmp-pharma.it email : info@cmp-pharma.it

5.0 Summary of modifications required if any

Item	Description	Action required	Projected completion date
1	Outfeed trays labels have are white.	Change labels color (GREEN for accepted and RED for rejection)	Before delivering
2		Cmp will send pictures of doors as in the layout	Before delivering
3	Too less space between outfeed Trays	Have outfeed trays space between trays (2cm)	Before delivering
4			
5			
6			
7			

6.0 Dates and Deadlines

	Planned Date	CMP signature	CUSTOMER signature
Delivery	to be defined	<i>David Tomer</i>	<i>张岳峰 2017.08.31</i>
Installation / commissioning	to be defined	<i>David Tomer</i>	<i>张岳峰 2017.08.31</i>

Comments

FAT RS1	TFDA	Serial No: 618	Pag. 16 - 17
---------	------	----------------	--------------

7.0 Approval

The FAT was conducted in the prescribed manner and to level of detail required and all actions or results/outcomes documented.



The Factory Acceptance Test is complete and accepted without reserves:

Name/Company	Signature	Date
CMP PHARMA	<i>Dante Poma</i>	2017-08-31
TFDA	<i>张磊</i>	2017.08.31
Name/Company	Signature	Date
Name/Company	Signature	Date


FAT RS1	TFDA	Serial No: 618	Pag. 17 - 17
---------	------	----------------	--------------











CMP PHAR.MA S.R.L. - via G. Galilei n° 26 36030 Costabissara (Vicenza) Italy - tel. 0039.0444.971716 - fax: 0039.0444.553262 -
 C.F.P.IVA 03972190247 BANCA POPOLARE DI VICENZA - SWIFT BPVI IT 21040 - IBAN IT 91Y0572860210040570004133
 website : www.cmp-pharma.it email : info@cmp-pharma.it

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