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Fibroblast culture from bone marrow

Materials:

- 10cm plates
- ChangD media (Irvine Scientific Cat. No. T105-500ml)
- Penicillin/Streptomycin (Gibco Cat. No. 15140-122) 用於未長出 colony 前 · 加進 ChangD
- Trypsin (0.05% EDTA)
- 1X PBS
- Freeze media (90% FBS + 10% DMSO)
- Grow in a hypoxic incubator for best results (5% O₂, 5% CO₂, 37 C). It can also be cultured in a regular normoxic incubator, but it will grow a bit slower

Procedure:

In ward: Collect 100-200µl BM in 5-10ml N/S (PBS) in sterile container (urine cup)

In lab

1. pipette above to 15ml conical tube, wash urine cup once with 3ml PBS.
2. centrifuge in RT (room temperature), 1200-1500rpm, 5min.
3. suction off supernatant, 可留一點點 N/S, 先用手敲打散 pellet · resuspend pellet with (ChangD + 1% Pen/Strep) 8-10ml, and plate on a 10cm tissue culture plate
4. 24-36 hrs later, **change media**: suction off media and replace with (ChangD + 1% Pen/Strep) 10ml
5. check colonies grow status QD or QOD, **change media** once to twice a week
6. if you have a fair amount of colonies (70-80% confluent), split cells: trypsinize and replat

Estimation

D4-6: individual spindle-shaped fibroblast

D7-10: small clusters


*wait until D14 for fibroblasts grow, if still no growth: failed

W3-5: 70-80% confluent

Splitting Primary Fibroblasts

When it's around 70-80% confluent, you can split cells

1. Aspirate media (plate 斜置) and wash cells with 5 mL of PBS (沿左側邊緣).
2. Aspirate PBS (plate 斜置、於右側抽吸), add 1-2 mL of trypsin (用滴的), and incubate in the cell incubator for at least 1 min (讓細胞游離).
3. Add each 9 ml of Chang D media to 2-6 new plates (分盤的數量取決於細胞長的速度)
4. Check that cells are off the plate (細胞已游離、呈圓形) and add 3-4 mL of Chang D media to plate. **Break cells apart** (反覆抽吸細胞、確定黏在 plate 的細胞都被洗下來).
5. Collect all cells and transfer to a 15ml conical tube. Spin down at 1500rpm for 5 minutes
6. Suction off supernatant, flick the tube to loosen up cell pellet and resuspend with 1 mL of Chang D media.
7. Do a 1:8 split if cells are growing really fast by adding 125 ul of the resuspended pellet to the 8 plates. If cells are growing slowly, do a 1:2 plate into both plates by adding 500 ul to both plates. Give cells a shake (前後左右、勿轉圈) and incubate.

Freeze fibroblasts

When it's around 70-80% confluent, you can freeze down to store.

Freeze at least **5 vials** each time.

1. Aspirate media (plate 斜置) and wash cells with 5 mL of PBS (沿左側邊緣).
2. Aspirate PBS (plate 斜置、於右側抽吸), add 1-2 mL of trypsin (用滴的), and incubate in the cell incubator for at least 1 min (讓細胞游離).
3. Check that cells are off the plate (細胞已游離、呈圓形) and add 3-4 mL of Chang D media to plate. **Break cells apart** (反覆抽吸細胞、確定黏在 plate 的細胞都被洗下來).
4. Collect all cells and transfer to a 15ml conical tube. Spin down at 1500rpm for 5 minutes
5. Suction off supernatant, flick the tube to loosen up cell pellet and resuspend with 1 mL (FBS + 10% DMSO). Aliquot in 1ml in cryotubes
6. Place cryotubes in cryo container and store at -80C overnight. Transfer tubes to liquid nitrogen storage

Label

Sample No.	AML-001
BM-F P0	BM-F P0
1x10cm plate	0.5x10cm plate
Processed date, Initials	2022/10/31 ML

Primary fibroblast culture from BM: P0 (passage 0)

Split with trypsin: P1, P2...

Generation of Fibroblasts from BM-MNC (Bone Marrow Mononuclear Cells)

Protocol:

1. Prepare before taking cells out of LN2 (liquid nitrogen): equilibrate Chang D media and 50 mL conical of FBS to RT (okay if cold). Take BM MNC cryovial from liquid nitrogen and place on dry ice.
2. Put the cryovial in bead or water bath and let thaw until there is only a small chunk of ice remaining. 幾乎融但未完全融化 · 進 hood 前先酒精消毒外面
3. Put 5mL of FBS into 15mL conical for resuspension/ centrifugation.
4. Using a 2mL serological pipette, take up 1mL FBS and gently resuspend cells in cryovial
 - a. **Gently and quickly** (max 2 tubes) repipette 2 or 3 times until the ice is fully melted.
5. Take up 2mL of cells from the cryovial and gently pipette dropwise on top of FBS in 15mL conical. 用滴的
6. Spin down for 5 mins at 1200 rpm at RT (Wash #1).
7. Pipette Chang D media 8-10 mL into labeled 10cm plate.
 - a. Label with Sample No., BM-F with passage number, thaw date, and initials.
8. When cells are done spinning, aspirate media.
 - a. Don't get too close to the pellet because it will be loosely packed.
9. Resuspend cells in 1 mL Chang D. Plate 300 uL of the resuspended cells on 10 cm plate to grow fibroblasts. **About remaining 700 uL:**
 - a. separate into 2-3 more plates for generating fibroblast
 - b. Add 1mL PBS to the remaining 700uL of BM MNC cells in the conical.
 - i. Count total cells and store back or
 - ii. Proceed to DNA extraction (somatic mutation)
10. Change Chang D media the following day then twice a week until ready (70-80% confluent) for split. Eventually, 1 plate to extract DNA and 5 to freeze back.

Fibroblasts growth from BM-F

Protocol: steps 1-8 is same as above

9. Resuspend cells in 1 mL Chang D. Plate all resuspended cells on 10 cm plate to grow fibroblasts.
10. Change Chang D media the following day then maybe ready (70-80% confluent) for split 3 days later (usually within 1 week).

DNA Extraction

[Different steps]

BM-MNC or BM-Fibroblast samples

DNeasy Blood & Tissue Kit (Ref # 69506) (Collection tube: DNeasy Mini Spin Column)

1. Label Epp tubes (1, 2, 3...) and have corresponding subject ID # recorded in lab notebook. **Note:** Maximum number of cells to load a column is 5×10^6 . In general, it is okay to process one 10cm plate of patient primary fibroblasts per column.
2. Resuspend (thawed) cell pellet in 200ul PBS.
3. Add 20ul **proteinase K** (in the kit box at RT).

Blood samples

QI Amp DNA Blood Mini Kit (Ref # 51106)

1. Label Epp tubes (1, 2, 3...) and have corresponding subject ID # recorded in lab notebook. Process 1 column per sample at a time.
2. Add 20ul **QIAGEN Protease** into the bottom of an Epp tube. Quiagen Protease kept in 4C.
3. Add 200ul blood per tube (can add PBS to bring volume up to 200ul) (keep remainder of blood in vial on ice since will only process 200ul at a time. Flash freeze vial after confirm extracted enough DNA. Otherwise, bring vial to RT again and process another 200ul aliquot).

[Same steps]

4. Add 200ul AL buffer, shake 15 seconds, then vortex for 30s, and quick spin in Mini Mouse.
5. Place Epp tubes in heatblock set to 56C for 10 minutes.
6. Spin down Epp tubes in Mini Mouse.
7. Add 200ul ETOH (100%, 200 proof), shake 15 seconds, vortex 30s to mix thoroughly, and quick spin.
8. Transfer to spin column in 2ml collection tube, close the cap.
9. Centrifuge at 6,000g for 1 min, then transfer column to a fresh collection tube.
10. Add 500ul buffer AW1, centrifuge at 6,000g for 1min, then transfer column to a fresh collection tube.
11. Add 500ul buffer AW2, centrifuge at 17,000g for 3 min, then transfer column to a fresh collection tube.
12. Dry spin, 17,000g for 1 min.
13. Transfer column to a clean Epp tube (labeled on the side with Sample ID # and today's date).
14. Add AE (**PB: 80ul; BM-F or BM-MNC: 100ul**) to the column, incubate for 5 min, then centrifuge at 6,000g for 1 min.
15. Pipette the eluted DNA back on the column to re-elute, incubate for 5 min, then centrifuge at 6,000g for 1 min.
16. Nanodrop (blanked with AE), record DNA concentration (ng/ul), 260/280, and 260/230, and total eluted volume.
17. Print labels for tubes, store in -20 freezer and document on the -20 C log.
Top = ID #. Side Labels:

Subject ID #

BL / BM-F / BM-MNC DNA

Date collected (collt'd ____)

DNA ng/ul

Date processed & initials


<Sanger sequence PCR primer design>

chr17:74737110A>C, RAB37, NM_001163989.1:c.101A>C (p.Asn34Thr), exon 1
 NC_000017.11 (74671131..74747335)

Ensembl

<https://useast.ensembl.org/index.html>

choose yellow one but **with interested position**

Choose exon

Configure Page | Personal Data

Display options

Select from available configurations: Current unsaved | Save current configuration

Display options

Flanking sequence at either end of transcript: 150

Number of base pairs per row: 60 bp

Intron base pairs to show at splice sites: 150

Show full intronic sequence:

Show exons only:

Line numbering: None

Show variants: Yes

Hide variants longer than 10bp:

Hide variants by frequency (MAF): Don't hide

Filter variants by consequence type:

- No filter
- 3 prime UTR variant
- 5 prime UTR variant
- NMD transcript variant

Filter variants by evidence status: No filter

Exons/ Introns | Translated sequence | Flanking sequence | Intron sequence | UTR

Variants

- 3 prime UTR
- 5 prime UTR
- Coding sequence
- Downstream
- Frameshift
- Inframe deletion
- Intronic
- Missense
- Splice acceptor
- Splice donor
- Splice region
- Start retained
- Stop gained
- Synonymous
- Upstream

Markup loaded

Sequence

```
.....tgagctcgggcgggtgacttacagatga
ggaagtgtctcggggccatcgggggagggcgagcccaggggtcccagctcccgcctcg
ccacctggggacagcccacggccgggctcggggcccgctgctgctcgggtgcgagc
GACTACGGGAACCTTCOCGAGCAGACGGGGTCCCCGCCGCCGCTCCCCAGGGGCAAG
CAAGCGACCACAGGGGACCGGTCCCGGGCTGGATGTGGCTCATGTCCGAAGCGCACGGA
GCCGAGCCGGTGTGCTCAGGGAGGCTGCCGCCCTTCACGCAGACCCCTGGGCTCTGC
GTGCCCTCAGGGAACAGCAAG
gtcggagcgggtctcgtcagggggcgacgggacggaggaggagcctgaggggtcccgg
tcggaggaggaggaggtggcggggggggtggggcgttctgagctctctctctcggc
tcggggccgggactgctcacctctcgtcca.....ggcggtggcttctggaagcaggtaggcct
ttgtgtcagcattcgtgccagccccctcttctcagatcctctccatgtgtctctct
ccttgaatcccagaagctcccctgaccccattactgctctgccccagcccctag
GTGATGCTTCTGGGAGACAGCGCGTGGGAAAACATGTTTCTCATCCAATTGAAGAG
GGGGCCTTCTCTCCGGAACTTCATAGCCACCGTCGGCATAGACTTCAGG
ctgaggtggtcagggacttgcttcagcagagagccagggtctggtcagggcagggg
gggttgccccaccttgctcacctggctcccagggactccagggctcagctcggag
ggcacaacecgctccccagaccacag.....ctgctgcccctgatggtatgatctggctgga
gactgtctgggctactgcacccactctaggcctggagagggaacaagaaggagcgtc
tcagagctgagagccacatgactcctgcccctccatccctggtcttttctcttccag
AACAGGTGGTGACTGTGGATGGGTGAGCTGAGCTGCAAG
gtgagaacagaggctggagttggggagggaaggatggaggacctgctccttccctccacc
tgaaccacaggaggtctgagccctgcctccgctggggcaattctgtgaggccacac
gggaggaatggcttttggttatttggadat.....atgtgatgctcaggcctcagaaacctt
aagaagaaaacgcccaccttcccacagatctcatacaacaagagcggcctgtccca
cccagcactctctcagcctgctcttccaccatttctccatcccaccctccag
```

● Double check primers

Primer BLAST

<https://www.ncbi.nlm.nih.gov/tools/primer-blast/>

The screenshot shows the Primer BLAST web interface with the following settings:

- Primer Parameters:**
 - Use my own forward primer (5'->3' on plus strand):
 - Use my own reverse primer (5'->3' on minus strand):
 - PCR product size: Min Max
 - # of primers to return:
 - Primer melting temperatures (T_m): Min Opt Max Max T_m difference
- Exon/intron selection:**
 - Exon junction span:
 - Exon junction match: Min 5' match Min 3' match Max 3' match
 - Intron inclusion: Primer pair must be separated by at least one intron on the corresponding genomic DNA
 - Intron length range: Min Max
- Primer Pair Specificity Checking Parameters:**
 - Specificity check: Enable search for primer pairs specific to the intended PCR template
 - Search mode:
 - Database: (highlighted in yellow)
 - Exclusion: Exclude predicted Refseq transcripts (accession with XM, XR prefix) Exclude uncultured/environmental sample sequences

Note: Parameter values that differ from the default are highlighted in yellow.

"Refseq representative genomes" that may contain these ALT_REF_LOCI records; "Genomes for selected eukaryotic organisms (primary assemblies only)" should not contain these. Including those records is a more comprehensive search, but since alignments to those records is almost always non-informative, you may decide to avoid them for convenience.

Primer pair 1

	Sequence (5'→3')	Length	Tm	GC%	Self complementarity	Self 3' complementarity
Forward primer	AAGTGTCTCGGGCCATC	18	59.01	61.11	4.00	1.00
Reverse primer	CCACTCCTCCCCTCCCTC	18	60.04	72.22	1.00	0.00

Products on target templates
 >NC_000017.11 Homo sapiens chromosome 17, GRCh38.p14 Primary Assembly

product length = 399
 Features associated with this product:
 ras-related protein rab-37 isoform 3
 ras-related protein rab-37 isoform 6

Forward primer 1 AAGTGTCTCGGGCCATC 18
 Template 74736799 74736816

Reverse primer 1 CCACTCCTCCCCTCCCTC 18
 Template 74737197 74737180

product length = 477
 Features flanking this product:
 218917 bp at 5' side: uncharacterized protein loc124904103
 222033 bp at 3' side: trinucleotide repeat-containing gene 6c protein isoform x13

Reverse primer 1 CCACTCCTCCCCTCCCTC 18
 Template 77782173 ..G....T..... 77782156

Reverse primer 1 CCACTCCTCCCCTCCCTC 18
 Template 77781697T.T..... 77781714

product length = 1222
 Features flanking this product:
 113792 bp at 5' side: protein sidekick-2 isoform x4
 444855 bp at 3' side: 60s ribosomal protein 138

- Double check primers

BLASTn: IDT protocol

<https://blast.ncbi.nlm.nih.gov/Blast.cgi>

<https://www.idtdna.com/pages/education/decoded/article/tips-for-using-blast-to-locate-pcr-primers>

Use these tips to refine Primer-BLAST results:

1. Concatenate the two primer sequences into one sequence separated by 5–10 Ns and enter into BLAST sequence box.
TGCTGTCTCCCATTTCCAGTNNNNNNNNNNCACTACCTCTACCCTCCGTT
2. Before submitting, narrow the search by selecting the species, if known; otherwise, choose **Nucleotide Collection (nr/nt)**. If you’re looking for RT-PCR primers, select the reference **mRNA sequences (refseq_mRNA)** database.
3. Under Program Selection, select the **Somewhat similar sequences (blastn) program**.
4. Under **Algorithm parameters**, decrease **word size to 7**, increase **expect threshold to 1000**, and **turn off the low complexity filter**.

In the example below, the results with a line connecting the 2 boxes indicate the 2 primers are in the same sequence (Figure 1A). Clicking on these results shows location within the sequence (between bases 1005–1766 of the *MRE11* transcript in yeast) and indicate the expected PCR product should be 762 bp in length (Figure 1B).


Figure 1. BLAST results show primer positioning and PCR fragment size. (A) Only the best alignment is shown (blue lines at the bottom). (B) Selecting the alignment provides the actual sequence with base by base alignment of the 2 portions of sequence showing 100% identity.

BLAST® » blastn suite Home Recent Results Saved Strategies Help

Take the BLAST survey today [Start survey](#)

Standard Nucleotide BLAST

blastn | blastp | blastx | tblastn | tblastx

BLASTN programs search nucleotide databases using a nucleotide query. more... [Reset page](#) [Bookmark](#)

Enter Query Sequence

Enter accession number(s), gi(s), or FASTA sequence(s) [Clear](#) Query subrange [?](#)

TGCTGTCTCCATTCCAGTNNNNCACTACCTCTACCCTCCGTT From

To

Or, upload file No file chosen [?](#)

Job Title [?](#)

Enter a descriptive title for your BLAST search [?](#)

Align two or more sequences [?](#)

Choose Search Set

Database Standard databases (nr etc.) rRNA/ITS databases Genomic + transcript databases Betacoronavirus

Experimental databases [Try experimental taxonomic nt databases](#) [Download](#)

For more info see [What are taxonomic nt databases?](#)

Organism [?](#)

Optional exclude [Add organism](#)

Enter organism common name, binomial, or tax id. Only 20 top taxa will be shown [?](#)

Exclude Models (XMP) Uncultured/environmental sample sequences

Optional Sequences from type material

Optional [YouTube](#) [Create custom database](#)

Optional Enter an Entrez query to limit search [?](#)

Program Selection

Optimize for Highly similar sequences (megablast)

More dissimilar sequences (discontinuous megablast)

Somewhat similar sequences (blastn)

Choose a BLAST algorithm [?](#)

[Feedback](#)

Algorithm parameters [Restore default search parameters](#)

General Parameters

Max target sequences [?](#)

Select the maximum number of aligned sequences to display [?](#)

Short queries Automatically adjust parameters for short input sequences [?](#)

Expect threshold [?](#)

Word size [?](#)

Max matches in a query range [?](#)

Scoring Parameters

Match/Mismatch Scores [?](#)

Gap Costs [?](#)

Filters and Masking

Filter Low complexity regions [?](#)

Species-specific repeats for: [?](#)

Mask Mask for lookup table only [?](#)

Mask lower case letters [?](#)

BLAST Search database Nucleotide collection (nr/nt) using Blastn (Optimize for somewhat similar sequences)

Show results in a new window

Download ▾ GenBank Graphics Sort by: E value ▾ Next Previous Descriptions

Homo sapiens chromosome 6 genomic scaffold, GRCh38.p14 alternate locus group ALT_REF_LOCI_7 HSCHR6_MHC_SSTO_CTG1
 Sequence ID: [NT_167249.2](#) Length: 4929269 Number of Matches: 6

Range 1: 2023626 to 2023645 [GenBank](#) [Graphics](#)  Next Match Previous Match

Score	Expect	Identities	Gaps	Strand
37.4 bits(40)	1.4	20/20(100%)	0/20(0%)	Plus/Plus

Query 1 TGCTGTCTCCCATTTCCAGT 20
 Sbjct 2023626 TGCTGTCTCCCATTTCCAGT 2023645

Range 2: 2023945 to 2023964 [GenBank](#) [Graphics](#) Next Match Previous Match First Match

Score	Expect	Identities	Gaps	Strand
37.4 bits(40)	1.4	20/20(100%)	0/20(0%)	Plus/Minus

Query 32 CACTACCTCTACCCTCCGTT 51
 Sbjct 2023964 CACTACCTCTACCCTCCGTT 2023945

Range 3: 4878141 to 4878154 [GenBank](#) [Graphics](#) Next Match Previous Match First Match

Score	Expect	Identities	Gaps	Strand
26.5 bits(28)	2610	14/14(100%)	0/14(0%)	Plus/Plus

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
NCBI Reference Sequence: NT_167249.2

[GenBank](#) [FASTA](#)

Read about these ALT_REF_LOCI records here:

<https://www.ncbi.nlm.nih.gov/grc/help/faq/#difference-between-alternate-loci-and-novel-patch>

Alternate loci and novel patches enable the reference assembly to represent allelic diversity. They are scaffold sequences that are given chromosome context through alignments to the corresponding chromosome regions. Alternate loci scaffolds and their alignments are included in major assembly releases, while novel patch scaffolds and their alignments are included in subsequent patch releases for that assembly. They can be considered functionally equivalent, as novel patches will be reassigned to the role of alternate loci scaffolds at the time of the next major assembly release. Assembly regions for which the GRC provides alternate loci or novel patch scaffolds are typically those with known alternate haplotypes (e.g. immune-associated regions), highly variable genomic regions (e.g. olfactory receptor regions) or those where there are structural variants having 5 Kb or more sequence not represented on the chromosome. Human alternate loci and all novel patch scaffolds also include one or more anchor sequence components to ensure their robust alignment to the chromosomes. Anchor sequences are component(s) that are also found in the corresponding chromosome. The sequence locations corresponding to anchor components are annotated on the GenBank records for all alternate loci and patch scaffolds. For more detail on patches please see [Introductions to Patches](#).