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


Figure 2

100µl virus diluent/well + 100µl cell suspension

5.4 Calculation of the virus titre by the Kärber formula

$\log \text{CCID}_{50} = L - d(S - 0.5)$, where:

L = log of lowest dilution used in the test

d = difference between log dilution steps

S = sum of proportion of "positive" tests (i.e. cultures showing CPE)

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
Poliovirus Isolation New Algorithm

- a Observed for a minimum of 5 days
- b Observe until $\geq 3+$ CPE obtained (usually 1-2 days, 5 days maximum; re-inoculate when toxicity or contamination observed)
- c Total minimum observation time of 10 days (2 x 5 days)
- d Pool positive tubes (if both tubes show $> 3+$ CPE on the same day) before final RD passage
- e Isolates can be serotyped by laboratories with an interest in NPEV diagnosis or to confirm proficiency

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型內鑑定(Intratypic differentiation, ITD)

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
RT-PCR condition (Do not use Ramp. The following in a normal cycler is best.

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Sequence Reporting

- The trimmed VP1 sequence is then reported using the following criteria:
 - Sequences classified as serotype 1 but having less than 85% identity (~136nt diff) with Sabin1 are designated as **PV1 wild***.
 - Sequences classified as serotype 1 and ≤ 9 nucleotide differences from Sabin1 are designated as **PV1 Sabin**.
 - Sequences classified as serotype 1 and ≥ 10 nucleotide differences from Sabin1 are designated as **PV1 VDPV**.
 - Sequences classified as serotype 2 and ≤ 5 nucleotide differences from Sabin2 are designated as **PV2 Sabin**.
 - Sequences classified as serotype 2 and ≥ 6 nucleotide differences** from Sabin1 are designated as **PV2 VDPV*****.
 - Sequences classified as serotype 3 but having less than 85% identity (~136nt diff) with Sabin3 are designated as **PV3 wild***.
 - Sequences classified as serotype 3 and ≤ 9 nucleotide differences from Sabin3 are designated as **PV3 Sabin**.
 - Sequences classified as serotype 3 and ≥ 10 nucleotide differences from Sabin3 are designated as **PV3 VDPV**.
- All chromatograms and consensus sequences should be archived and catalogued in a manner to facilitate easy retrieval of raw sequence data and analysis of completed consensus data.

**Exceptions to this rule may occur, including the possibility of a VDPV that has replicated for more than 10-15 years. Additional analyses can be performed to determine if a sequence has a distant Sabin link.*

*** Determined at the Annual GPLN September 2010 Meeting for PV2 VDPV definition*

**** Please note that PV2 wild virus is not listed because no indigenous PV2 wild virus has circulated since 1999. It is important to remember that a wild PV2 virus might be detected if cross-contamination by type 2 reference strains has occurred.*

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